




LAB REPORT

Abrasion Test Results

The abrasion test employed two bar samples measuring 1" x 3" x 1/4". One of the samples was Kryptane 8008 and the other samples were various materials to be tested against Kryptane. Both materials were rotated at 3900 RPM for 20 hours in a 50/50 - 20 grit silica sand and water slurry at 100 degrees F. Volume loss from each of the eight samples were computed and compared. Included in the test were the most popular abrasion resistant materials available. The other urethanes wore more than three times faster than Kryptane 8008!

<u>SAMPLE</u>	<u>RELATIVE VOLUME LOSS</u>
KRYPTANE 8008	10
Urethane #1	48
Urethane #2	32
Urethane #3	33
Urethane #4	174
Urethane #5	61
Hardened Tool Steel	270
UHMW Polyethylene	112
Natural Rubber	58



The information in this bulletin is derived from the best available sources and is believed to be accurate. However, no guarantee is expressed or implied regarding the accuracy of this data or the use of this product, nor are any statements in this bulletin intended to infringe on any patent.